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Measurement and instrumentation

**Ref :
281 262**

Français – p 1

English – p 5

Version : 6009

Alimentation stabilisée EL5
EL5 stabilised power supply

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
- Power supply protected against all operating errors.
- Increased safety due to double insulation.
- Suitable for use in classes and for experiments.
- Electronic protection.
- Stabilised DC voltage.

1 Purpose


EL5 is a stabilised power supply designed for electricity and electronic experiments.

EL5 is connected to the 220/230 volts mains and supplies the user with:

- A 6 or 12 volts DC voltage up to 1.5 amperes. Electronic protection.
- A 6 or 12 volts AC voltage up to 5 amperes. Protection: thermal bi-metal switch.

2 Commissioning

1. Connect EL5 to the 220/230 V mains supply.
2. Place the On/Off switch on \perp . The indicator lamp comes on.
3. Select the desired voltage with the 6 - 12 V selector

4. Select DC or AC with the  or  selector
5. Connect your circuit (lamp, motor, etc.) on the DC or AC sockets.

Notes:

- In DC mode, always connect the circuit after powering on the EL5. A sudden surge of current can actually trigger the safety device (i instantaneous > 1.5 A) and EL5 may appear not to function. In this case the circuit must be disconnected and reset.
- The electronic circuit of the power supply generates heat, so do not block the ventilation vents.

3 Technical specifications

- Mains power supply 230 V, $\pm 5\%$.
- DC 6 - 12 V / 1.5 A :
 - Voltage drop < 50 mV for a load between 0 and 1.5 A.
 - Residual ripple < 50 mV peak to peak.
- AC 6 - 12 V / 5 A :
 - Voltage drop < 10 % for a load between 0 and 5 A.
 - Dimensions: 150 x 110 x 170 mm.
 - Mass: 2.9 kg.
 - Primary protection: Non accessible tubular fuse
 - AC output protection: thermal circuit breaker with quick cut off.
 - Auto-reset.
 - DC output protection: electronic limitation. Auto-reset after cut off of the user circuit (disconnect a wire from the circuit or insert a switch in the circuit).
- This device complies with EN 61010-1, installation category II, pollution degree 2.

4 Care – maintenance – repairs

4.1 Systematic research in the case of a malfunction

If there is no output (zero current) after connecting the circuit to be used to the power supply:

1. The output voltage is zero:
 - Check that the selector switch located between AC and DC is properly positioned.
2. The output voltage is not zero but drops to zero after connecting the circuit:
 - Check your loading circuit (maybe the current in the circuit exceeds 5 A, or there is a short circuit, etc.).

In AC, the power supply is protected by a thermal circuit breaker, which resets automatically with a small click.


- Disconnect your circuit and make the necessary correction: connection to be corrected, rheostat to be placed in series, etc.

5 After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service * :
+ 33 (0)2 32 29 40 50

** Call cost not included*